DOC316.53.01109

Phosphonates

PhosVer 3 with Persulfate UV Oxidation¹ Multiple ranges from 0.02 to 125 mg/L PO₄³⁻

Method 8007 Powder Pillows

Scope and application: For boiler and cooling water, wastewater and seawater.

Adapted from Blystone, P., Larson, P., A Rapid Method for Analysis of Phosphate Compounds, International Water Conference, Pittsburgh, PA. (Oct 26-28, 1981)



Test preparation

Instrument-specific information

Table 1 shows all of the instruments that have the program for this test. The table also shows sample cell and orientation requirements for reagent addition tests, such as powder pillow or bulk reagent tests.

To use the table, select an instrument, then read across to find the applicable information for this test.

Table 1 Instrument-specific information

| Instrument | Sample cell orientation | Sample cell |
|------------|--|-------------|
| DR 6000 | The fill line is to the right. | 2495402 |
| DR 3800 | | |
| DR 2800 | | 10 mL |
| DR 2700 | | |
| DR 1900 | | |
| DR 5000 | The fill line is toward the user. | |
| DR 3900 | | |
| DR 900 | The orientation mark is toward the user. | 2401906 |

Before starting

Install the instrument cap on the DR 900 cell holder before ZERO or READ is pushed.

Clean all glassware with 6.0 N (1:1) hydrochloric acid, then fully rinse with deionized water to remove contaminants.

Do not use a detergent that contains phosphate to clean glassware. The phosphate in the detergent will contaminate the sample.

Wear UV safety goggles while the UV lamp is on.

Do not touch the UV lamp surface with bare fingers. Fingerprints can damage the glass. Rinse the lamp and wipe with a soft, clean tissue between tests.

The UV digestion in this procedure is normally complete in less than 10 minutes. However, high-organic loaded samples or a weak lamp can cause incomplete phosphate conversion. To check conversion efficiency, use a longer digestion time and make sure the readings do not increase.

Review the Safety Data Sheets (MSDS/SDS) for the chemicals that are used. Use the recommended personal protective equipment.

Dispose of reacted solutions according to local, state and federal regulations. Refer to the Safety Data Sheets for disposal information for unused reagents. Refer to the environmental, health and safety staff for your facility and/or local regulatory agencies for further disposal information.

Items to collect


| Description | Quantity |
|---|----------|
| Bottle, square, with 25-mL mark | 1 |
| Cylinder, graduated mixing, 50-mL | 1 |
| Goggles, UV safety | 1 |
| Pipet, serological, 10-mL | 1 |
| PhosVer® 3 Phosphate Reagent Powder Pillows, 10-mL | 2 |
| Potassium Persulfate Powder Pillow for Phosphonate | 1 |
| Pipet filler, safety bulb | 1 |
| Sample cells (For information about sample cells, adapters or light shields, refer to Instrument-specific information on page 1.) | 2 |
| Water, deionized | varies |
| UV lamp with power supply | 1 |

Refer to Consumables and replacement items on page 7 for order information.

Sample collection and storage

- Collect samples in clean glass or plastic bottles that have been cleaned with 6 N (50%) hydrochloric acid and rinsed with deionized water.
- Do not use a commerical detergent to clean the sample bottles. The phosphate in the detergent will contaminate the sample.
- To preserve samples for later analysis, adjust the sample pH to less than 2 with concentrated sulfuric acid (approximately 2 mL per liter). No acid addition is necessary if the sample is tested immediately.
- Keep the preserved samples at or below 6 °C (43 °F) for a maximum of 24 hours.
- Let the sample temperature increase to room temperature before analysis.
- Before analysis, adjust the pH to 7 with 5 N sodium hydroxide solution.
- Correct the test result for the dilution caused by the volume additions.

Powder pillow procedure with UV photolysis




1. Start program 501
Phosphonates. For
information about sample
cells, adapters or light
shields, refer to Instrumentspecific information
on page 1.

Note: Although the program name can be different between instruments, the program number does not change.



2. Select the sample size from Table 2 on page 5. Use a pipet to add the correct volume of sample into a 50-mL graduated cylinder.


If necessary, dilute the sample to 50 mL with deionized water and mix well.




3. Prepare the blank: Fill a sample cell to the 10-mL mark with the diluted sample from step 2.




4. Prepare the digested sample: Fill a mixing bottle to the **25-mL mark** with the diluted sample from step 2.



5. Add the contents of one Potassium Persulfate for Phosphonate Powder Pillow to the 25-mL sample.



6. Swirl to mix.



7. Put on UV safety goggles.



8. Put the ultraviolet lamp into the mixing bottle. Turn on the UV lamp.




9. Start the instrument timer. A 10-minute reaction time starts.

Phosphonates are converted to orthophosphate in this step.




10. When the timer expires, turn off the UV lamp. Remove the UV lamp from the sample.



11. Prepare the sample: Fill a second sample cell to the 10-mL mark with the digested sample.



12. Add the contents of one PhosVer 3 Phosphate Reagent Powder Pillow to both the blank and the prepared sample.




13. Immediately swirl both cells vigorously for 20–30 seconds to mix. Some powder may not dissolve.

A blue color shows if phosphate is present. Both the sample and the blank may show color.



14. Start the instrument timer. A 2-minute reaction time starts.


If the sample temperature is less than 15 °C (59 °F), wait 4 minutes for color development.




15. When the timer expires, clean the blank sample cell. Complete the rest of the steps in this procedure within 3 minutes.



16. Insert the blank into the cell holder.



17. Push **ZERO**. The display shows 0.00 mg/L PO_4^{3-} .



18. Clean the prepared sample cell.



19. Insert the prepared sample into the cell holder.



20. Push **READ**. Results show in $mg/L PO_4^{3-}$.



21. Multiply the results by the applicable sample volume multiplier in Table 2 on page 5 for the phosphonate concentration. Refer to Table 3 on page 5 to report results as the phosphonate compound.

Select the sample volume and multiplier

Use the expected phosphonate concentration to select a sample volume (refer to Table 2). Use the multiplier to adjust the test result (in mg/L PO_4^{3-}) for the sample volume that was used.

Table 2 Expected phosphonate range with multiplier

| Expected range (mg/L phosphonate) | Sample volume (mL) | Multiplier |
|-----------------------------------|--------------------|------------|
| 0–2.5 | 50 | 0.1 |
| 0–5 | 25 | 0.2 |
| 0–12.5 | 10 | 0.5 |
| 0–25 | 5 | 1 |
| 0–125 | 1 | 5 |

Convert phosphate to phosphonate

To convert the final test result from $mg/L PO_4^{3-}$ to active phosphonate, multiply the final test result by the applicable conversion factor in Table 3.

Table 3 Conversion factors by phosphonate type

| Phosphonate type | Conversion factor |
|------------------|-------------------|
| DETPMPA | 1.207 |
| EDTMPA | 1.148 |
| HEDPA | 1.085 |
| HMDTMPA | 1.295 |
| HPA | 1.49 |
| NTP | 1.05 |
| PBTC | 2.84 |

Interferences

Interference levels decrease as the sample size increases. For example, copper does not interfere at or below 100 mg/L for a 5.00 mL sample. If the sample volume is increased to 10 mL, copper will begin to interfere above 50 mg/L.

| Interfering substance | Interference level (5-mL sample) |
|--------------------------|---|
| Aluminum | 100 mg/L |
| Arsenate | Interferes at all levels |
| Benzotriazole | 10 mg/L |
| Bicarbonate | 1000 mg/L |
| Bromide | 100 mg/L |
| Calcium | 5000 mg/L |
| CDTA | 100 mg/L |
| Chloride | 5000 mg/L |
| Chromate | 100 mg/L |
| Copper | 100 mg/L |
| Cyanide | 100 mg/L (Increase the UV digestion to 30 minutes.) |
| Diethanoldithiocarbamate | 50 mg/L |
| EDTA | 100 mg/L |

| Interfering substance | Interference level (5-mL sample) |
|--|---|
| Iron | 200 mg/L |
| Nitrate | 200 mg/L |
| NTA | 250 mg/L |
| Orthophosphate | 15 mg/L |
| Phosphites and organophosphorus compounds | Reacts quantitatively. Metaphosphates and polyphosphates do not interfere. |
| Silica | 500 mg/L |
| Silicate | 100 mg/L |
| Sulfate | 2000 mg/L |
| Sulfide | Interferes at all levels |
| Sulfite | 100 mg/L |
| Thiourea | 10 mg/L |
| Highly buffered samples or extreme sample pH | Can prevent the correct pH adjustment of the sample by the reagents. Sample pre-treatment may be necessary. |

Accuracy check

Digestion method

To validate the full procedure with the digestion, prepare a solution that has a known concentration of a phosphonate compound. Use the test procedure to measure the concentration of the phosphonate solution.

Standard solution method

To validate the colorimetric portion of the procedure (without digestion), use a phosphate standard solution for the sample and deionized water for the blank. Add the PhosVer 3 reagent directly to 10 mL of the phosphate standard solution and to the blank. The expected result is 10 times the value of the standard solution due to a built-in dilution factor of 10 in the calibration.

Items to collect:

- Phosphate Standard Solution, 1 mg/L (the expected result is 10 mg/L if 10 mL is used)
- 1. Use the test procedure to measure the concentration of the standard solution.
- **2.** Compare the expected result to the actual result.

Note: The factory calibration can be adjusted slightly with the standard adjust option so that the instrument shows the expected value of the standard solution. The adjusted calibration is then used for all test results. This adjustment can increase the test accuracy when there are slight variations in the reagents or instruments.

Method performance

The method performance data that follows was derived from laboratory tests that were measured on a spectrophotometer during ideal test conditions. Users can get different results under different test conditions.

| Program | Standard | Precision (95% confidence interval) | Sensitivity Concentration change per 0.010 Abs change |
|---------|---|--|--|
| 501 | 2.00 mg/L PO ₄ ³⁻ | 1.97–2.03 mg/L PO ₄ ^{3–} | Refer to Sensitivity on page 7. |

Sensitivity

The sensitivity depends on the sample volume. Sensitivity is expressed as PO_4^{3-} in Table 4. To express as a specific phosphonate, refer to Table 3 on page 5.

Table 4 Sensitivity per sample volume

| Range (mg/L phosphonate) | Sample volume (mL) | Concentration change per 0.010 Abs change |
|--------------------------|--------------------|---|
| 0–2.5 | 50 | 0.02 mg/L PO ₄ ³⁻ |
| 0–5 | 25 | 0.04 mg/L PO ₄ ³⁻ |
| 0–12.5 | 10 | 0.10 mg/L PO ₄ ³⁻ |
| 0–25 | 5 | 0.20 mg/L PO ₄ ³⁻ |
| 0–125 | 1 | 1.00 mg/L PO ₄ ³⁻ |

Summary of method

This method is directly applicable to boiler and cooling tower samples. The procedure is based on a UV-catalyzed oxidation of phosphonate to orthophosphate. The orthophosphate reacts with the molybdate in the PhosVer 3 reagent to form a mixed phosphate/molybdate complex. This complex is reduced by the ascorbic acid in the PhosVer 3, which gives a blue color that is proportional to the amount of phosphonate in the original sample. The orthophosphate in the original sample is removed when the blank is used to set the zero concentration. The measurement wavelength is 880 nm (DR 1900: 710 nm) for spectrophotometers or 610 nm for colorimeters.

Consumables and replacement items

Required reagents

| Description | Quantity/test | Unit | Item no. |
|---|---------------|-----------|----------|
| Water, deionized | varies | 4 L | 27256 |
| Phosphonate Reagent Set, 10-mL | 1 | 100 tests | 2429700 |
| Includes: | | | |
| PhosVer® 3 Phosphate Reagent Powder Pillow ¹ , 10-mL | 1 | 100/pkg | 2106069 |
| Potassium Persulfate Powder Pillow for Phosphonate | 1 | 100/pkg | 2084769 |

¹ PhosVer is a registered trademark of Hach Company.

Required apparatus

| Description | Quantity/test | Unit | Item no. |
|---|---------------|------|----------|
| Bottle, square, with 25-mL mark | 1 | each | 1704200 |
| Beaker, polypropylene, 50-mL, low form | 1 | each | 108041 |
| Mixing cylinder, graduated, 50-mL, with glass stopper | 1 | each | 189641 |
| UV safety goggles | 1 | each | 2113400 |
| Pipet, serological, graduated, 10-mL | 1 | each | 53238 |
| Pipet filler, safety bulb | 1 | each | 1465100 |
| UV lamp with power supply, 115 VAC | 1 | each | 2082800 |
| OR | | | |
| UV lamp with power supply, 230 VAC | 1 | each | 2082802 |

Recommended standards

| Description | Unit | Item no. |
|--|--------|----------|
| Phosphate Standard Solution, 1-mg/L as PO ₄ ³⁻ | 500 mL | 256949 |

Optional reagents and apparatus

| Description | Unit | Item no. |
|---|----------|----------|
| Hydrochloric Acid Solution, 6 N (1:1) | 500 mL | 88449 |
| Sulfuric Acid, concentrated, ACS | 500 mL | 97949 |
| Thermometer, non-mercury, –10 to +225 °C | each | 2635700 |
| Paper, pH, 0–14 pH range | 100/pkg | 2601300 |
| Ampule Breaker, 10-mL Voluette [®] Ampules | each | 2196800 |
| PhosVer 3 Phosphate Reagent Powder Pillows, 10-mL | 1000/pkg | 2106028 |
| UV Lamp, shortwave, pencil type | each | 2671000 |
| Power supply, 115 V/60 Hz | each | 2670700 |
| Power supply, 220 V/50 Hz | each | 2670702 |
| Phosphate Standard Solution, 3-mg/L as PO ₄ ³⁻ | 946 mL | 2059716 |
| Phosphate Standard Solution, 10-mg/L as PO ₄ ³⁻ | 946 mL | 1420416 |
| Phosphate Standard Solution, 15-mg/L as PO ₄ ³⁻ | 100 mL | 1424342 |
| Phosphate Standard Solution, 30-mg/L as PO ₄ ³⁻ | 946 mL | 1436716 |
| Phosphate Standard Solution, 50-mg/L, 10-mL Voluette [®] Ampules | 16/pkg | 17110 |
| Phosphate Standard Solution, 100-mg/L as PO ₄ | 100 mL | 1436832 |
| Phosphate Standard Solution, 10-mL ampule, 500 mg/L as PO ₄ | 16/pkg | 1424210 |
| Phosphate Standard Solution, 500-mg/L as PO ₄ | 100 mL | 1424232 |