## DOC062.53.20128.Apr17

## **Parameter Specific Sample Preparation**

Parameters	COD/BOD <sub>5</sub>	тос	Inorganic nitrogen		N total	Phosphorus	Heavy metals	Chloride	Culmbata
			NH <sub>4</sub> -N NO <sub>2</sub> -N	NO <sub>3</sub> -N	N total	P <sub>total</sub>	total	Chloride	Sulphate
Short-term conservation	+2 to +5 °C (up to 24 hours)	+2 to +5 °C, pH <2 with sulphuric acid (up to 1 week)	Not recommended	+2 to +5 °C, pH <2 with sulphuric acid (up to 24 hours)	pH 3 to 7 (up to 1 week)	Not required	Not required	Not required	+2 to +5 °C (up to 1 week)
Longer-term conservation	-18°C (up to 14 days)	-18 °C, pH <2 (up to 1 month)	Not recommended		-18 °C (up to 14 days)	pH <2 with sulphuric acid (up to 1 month)	pH <2 with nitric acid (up to 1 month)	Up to 1 month	Not recommended
Preparation of sample	Homogenization (DIN 38 402, part A 30*)	Homogenization (DIN 38 402, part A 30*/33	Filtration if required		Homogenization (DIN 38 402, part A 30*)	Homogenization (DIN 38 402, part A 30*)	Homogenization (DIN 38 402, part A 30*/31), digestion	Filtration if required	Filtration if required
Sampling vessel material	Glass or polyethylene (PE)			Glass or polyethylene (PE) Borosilicate glass			Glass or polyethylene (PE)		

<sup>\*</sup>Homogenization of the sample for all parameters that are determined in the original sample (overall sample).

